

DICER1 Polyclonal Antibody



Catalog Number: 162826

Product name

DICER1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to DICER1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human DICER1 (NP_001258211.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - DICER1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using DICER1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.