

Defb26 Polyclonal Antibody



Catalog Number: 162891

Product name

Defb26 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to Defb26

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 22-128 of rat Defb26 (NP_001032595.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - Defb26 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using Defb26 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.