

C9orf72 Polyclonal Antibody



Catalog Number: 162941

Product name

C9orf72 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to C9orf72

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human C9orf72 (NP_060795.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Western blot - C9orf72 Polyclonal Antibody

Western blot analysis of extracts of mouse kidney, using C9orf72 antibody at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.