

# HNRNPCL1 Polyclonal Antibody



Catalog Number: 162972

---

## Product name

HNRNPCL1 Polyclonal Antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

IF

## Antibody description

Polyclonal antibody to HNRNPCL1


## Dilutions

IF 1:50 - 1:200

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human HNRNPCL1 (NP\_001013653.1).

## Validations



## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Immunofluorescence - HNRNPCL1 Polyclonal Antibody

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

Immunofluorescence analysis of U2OS cells using HNRNPCL1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.