

SLC18A3 Polyclonal Antibody



Catalog Number: 163000

Product name

SLC18A3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SLC18A3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SLC18A3 (NP_003046.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Immunohistochemistry - SLC18A3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using SLC18A3 antibody at dilution of 1:200 (40x lens).