

# IL23R Polyclonal Antibody



Catalog Number: 163041

## Product name

IL23R Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to IL23R

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 405-629 of human IL23R (NP\_653302.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Western blot - IL23R Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using IL23R antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.