

DTNBP1 Polyclonal Antibody



Catalog Number: 163143

Product name

DTNBP1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DTNBP1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human DTNBP1 (NP_115498.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Western blot - DTNBP1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DTNBP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - DTNBP1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using DTNBP1 antibody. Blue: DAPI for nuclear staining.