

ETV6 Polyclonal Antibody



Catalog Number: 163164

Product name

ETV6 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ETV6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ETV6 (NP_001978.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IP


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IP 1:50 - 1:100

Validations



Western blot - ETV6 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ETV6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunoprecipitation - ETV6 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of SW620 cells using 3ug ETV6 antibody . Western blot was performed from the immunoprecipitate using ETV6 antibody at a dilution of 1:500.