

# NUBP1 Polyclonal Antibody



Catalog Number: 163179

## Product name

NUBP1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to NUBP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 201-320 of human NUBP1 (NP\_002475.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Western blot - NUBP1 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using NUBP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - NUBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using NUBP1 antibody at dilution of 1:100



(40x lens).

# NUBP1 Polyclonal Antibody




Catalog Number: 163179

---


## Immunohistochemistry - NUBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human thyroid cancer using NUBP1 antibody at dilution of 1:100 (40x lens).




## Immunohistochemistry - NUBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using NUBP1 antibody at dilution of 1:100 (40x lens).



## Immunohistochemistry - NUBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using NUBP1 antibody at dilution of



1:100 (40x lens).

## Immunohistochemistry - NUBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using NUBP1 antibody at dilution of 1:100 (40x lens).