

SMN2 Polyclonal Antibody



Catalog Number: 163232

Product name

SMN2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SMN2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human SMN2 (NP_059107.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

IP 1:50 - 1:100

Validations



Immunoprecipitation - SMN2 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of 293T cells, using 3 ug SMN2 antibody . Western blot was performed from the immunoprecipitate using SMN2 antibody at a dilution of 1:1000.