

STK26 Polyclonal Antibody



Catalog Number: 163243

Product name

STK26 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to STK26

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 297-416 of human STK26 (NP_057626.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions


WB 1:500 - 1:2000

Validations



Western blot - STK26 Polyclonal Antibody

Western blot analysis of extracts of rat thymus, using STK26 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.



Western blot - STK26 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using STK26 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.