

ATP5I Polyclonal Antibody



Catalog Number: 163347

Product name

ATP5I Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ATP5I

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human ATP5I (NP_009031.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - ATP5I Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using ATP5I antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ATP5I Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using ATP5I antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - ATP5I Polyclonal Antibody

ATP5I Polyclonal Antibody

Catalog Number: 163347

Immunohistochemistry of paraffin-embedded human stomach using ATP5I antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - ATP5I Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using ATP5I antibody at dilution of 1:100 (40x lens).