

NR5A2 Polyclonal Antibody



Catalog Number: 163358

Product name

NR5A2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NR5A2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human NR5A2 (NP_003813.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:200

Validations



Immunohistochemistry - NR5A2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using NR5A2 antibody at dilution of 1:100 (40x lens).