

HOXD9 Polyclonal Antibody



Catalog Number: 163364

Product name

HOXD9 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to HOXD9

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HOXD9 (NP_055028.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions


WB 1:500 - 1:2000

Validations



Western blot - HOXD9 Polyclonal Antibody

Western blot analysis of extracts of SH-SY5Y cells, using HOXD9 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 5s.

Western blot - HOXD9 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HOXD9 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.