

# THBS1 Polyclonal Antibody



Catalog Number: 163392

## Product name

THBS1 Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to THBS1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human THBS1 (NP\_003237.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations




Immunohistochemistry - THBS1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse spinal cord using THBS1 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - THBS1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using THBS1 antibody at dilution of



1:200 (40x lens).

Immunohistochemistry - THBS1 Polyclonal Antibody

# THBS1 Polyclonal Antibody



Catalog Number: 163392

---

Immunohistochemistry of paraffin-embedded

mouse spleen using THBS1 antibody at dilution of 1:200 (40x lens).