

# UCP3 Polyclonal Antibody



Catalog Number: 163394

## Product name

UCP3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to UCP3

## Preparation

Antigen: Recombinant protein of human UCP3.

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Immunofluorescence - UCP3 Polyclonal Antibody


Immunofluorescence analysis of C6 cells using UCP3 antibody at dilution of 1:100. Blue: DAPI for



nuclear staining.

Immunofluorescence - UCP3 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using UCP3 antibody at dilution of 1:100. Blue: DAPI for



nuclear staining.

Immunofluorescence - UCP3 Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using UCP3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.