

HNRNPK Polyclonal Antibody



Catalog Number: 163398

Product name

HNRNPK Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to HNRNPK

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HNRNPK (NP_112552.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - HNRNPK Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HNRNPK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - HNRNPK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using HNRNPK Antibody at dilution of 1:100



(40x lens).



HNRNPK Polyclonal Antibody



Catalog Number: 163398

Immunohistochemistry - HNRNPK Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human esophagus using HNRNPK Antibody at dilution of 1:100 (40x lens).



(40x lens).

Immunohistochemistry - HNRNPK Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat liver using HNRNPK Antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - HNRNPK Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse heart using HNRNPK Antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - HNRNPK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using HNRNPK Antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - HNRNPK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using HNRNPK Antibody at dilution of 1:200

Immunohistochemistry - HNRNPK Polyclonal Antibody

HNRNPK Polyclonal Antibody



Catalog Number: 163398

Immunohistochemistry of paraffin-embedded

human esophagus using HNRNPK Antibody at dilution of 1:200 (40x lens).