

ARHGEF5 Polyclonal Antibody



Catalog Number: 163400

Product name

ARHGEF5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ARHGEF5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1248-1597 of human ARHGEF5 (NP_005426.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - ARHGEF5 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using ARHGEF5 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ARHGEF5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using ARHGEF5 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - ARHGEF5 Polyclonal Antibody

ARHGEF5 Polyclonal Antibody



Catalog Number: 163400

Immunohistochemistry of paraffin-embedded human stomach using ARHGEF5 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ARHGEF5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using ARHGEF5 antibody at dilution of 1:100 (40x lens).