

CADPS Polyclonal Antibody



Catalog Number: 163402

Product name

CADPS Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CADPS

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1100 of human CADPS (NP_899630.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - CADPS Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using CADPS antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CADPS Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CADPS antibody at dilution of 1:100 (40x lens).