

# ARL6IP4 Polyclonal Antibody



Catalog Number: 163435

## Product name

ARL6IP4 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Specificity

Human

## Antibody description

Polyclonal antibody to ARL6IP4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 313-402 of human ARL6IP4 (NP\_057722.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC


## Dilutions

## Validations



Immunohistochemistry - ARL6IP4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung using ARL6IP4 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ARL6IP4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using ARL6IP4 antibody at dilution of 1:100 (40x lens).