

# ASIC4 Polyclonal Antibody



Catalog Number: 163438

---

## Product name

ASIC4 Polyclonal Antibody

## Specificity

Rat

## Antibody description

Polyclonal antibody to ASIC4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 350-460 of human ASIC4 (NP\_878267.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions

IHC 1:50 - 1:200

## Validations



Immunohistochemistry - ASIC4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using ASIC4 antibody at dilution of 1:100 (40x lens).