

# TLR5 Polyclonal Antibody



Catalog Number: 163449

## Product name

TLR5 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TLR5

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 714-858 of human TLR5 (NP\_003259.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - TLR5 Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using TLR5 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.