

C12orf10 Polyclonal Antibody



Catalog Number: 163450

Product name

C12orf10 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to C12orf10

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 48-376 of human C12orf10 (NP_067653.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - C12orf10 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using C12orf10 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - C12orf10 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat ovary using C12orf10 antibody at dilution of 1:100



(40x lens).



Immunohistochemistry - C12orf10 Polyclonal Antibody

C12orf10 Polyclonal Antibody



Catalog Number: 163450


Immunohistochemistry of paraffin-embedded human lung cancer using C12orf10 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - C12orf10 Polyclonal Antibody

Immunohistochemistry - C12orf10 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using C12orf10 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human uterine cancer using C12orf10 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - C12orf10 Polyclonal Antibody

Immunohistochemistry - C12orf10 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using C12orf10 antibody at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded mouse testis using C12orf10 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - C12orf10 Polyclonal Antibody

C12orf10 Polyclonal Antibody



Catalog Number: 163450

Immunohistochemistry of paraffin-embedded

mouse kidney using C12orf10 antibody at dilution of 1:100 (40x lens).