

BRSK1 Polyclonal Antibody



Catalog Number: 163457

Product name

BRSK1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BRSK1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 310-510 of human BRSK1 (NP_115806.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - BRSK1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat liver using BRSK1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using BRSK1 antibody at dilution of 1:100




(40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

BRSK1 Polyclonal Antibody

Catalog Number: 163457


Immunohistochemistry of paraffin-embedded human stomach using BRSK1 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using BRSK1 antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded


mouse testis using BRSK1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using BRSK1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - BRSK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using BRSK1 antibody at dilution of 1:100 (40x lens).