

# Kctd19 Polyclonal Antibody



Catalog Number: 163474

---

## Product name

Kctd19 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to Kctd19

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse Kctd19 (NP\_001288102.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - Kctd19 Polyclonal Antibody

Immunofluorescence analysis of mouse testis using Kctd19 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.