

C10orf62 Polyclonal Antibody



Catalog Number: 163478

Product name

C10orf62 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to C10orf62

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human C10orf62 (NP_001009997.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - C10orf62 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using C10orf62 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - C10orf62 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using C10orf62 antibody at dilution of 1:100 (40x lens).