

Dimethyl-Lysine pAb



Catalog Number: 163490

Product name

Rabbit IgG

Dimethyl-Lysine pAb

Applications

Specificity

WB, IP

ALL

Dilutions

Antibody description

WB 1:500 - 1:1000

Polyclonal antibody to Dimethyl-Lysine


IP 1:20 - 1:50

Preparation

Validations

Antigen: Dimethyl Protein.

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.2.

Storage

- Dimethyl-Lysine pAb

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Dot-blot analysis of all sorts of methylation peptides using Dimethyl -Lysine antibody .

Polyclonal

Ig Type