

INHA Polyclonal Antibody



Catalog Number: 163495

Product name

INHA Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to INHA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-366 of human INHA (NP_002182.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - INHA Polyclonal Antibody

Western blot analysis of extracts of rat testis, using INHA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - INHA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using INHA antibody at dilution of 1:200 (40x lens).