

MYOG Polyclonal Antibody



Catalog Number: 163515

Product name

MYOG Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to MYOG

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-1-224 of human MYOG (NP_002470.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - MYOG Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using MYOG antibody at dilution of 1:100. Blue: DAPI for nuclear staining.