

# TFF1 Polyclonal Antibody



Catalog Number: 163552

## Product name

TFF1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TFF1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-84 of human TFF1 (NP\_003216.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:200

## Validations




## Western blot - TFF1 Polyclonal Antibody

Western blot analysis of extracts of BT-474 cells, using TFF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



## Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human mammary cancer using TFF1 antibody at dilution of 1:200 (40x lens).




# TFF1 Polyclonal Antibody

Catalog Number: 163552


## Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TFF1 antibody at dilution of 1:200 (40x lens).



## Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using TFF1 antibody at dilution of



1:200 (40x lens).

## Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat lung using TFF1 antibody at dilution of 1:200 (40x



lens).


## Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TFF1 antibody at dilution of 1:200 (40x lens).



## Immunohistochemistry - TFF1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using TFF1 antibody at dilution of 1:200



(40x lens).

## Immunohistochemistry - TFF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TFF1 antibody at dilution of 1:200 (40x lens).



# TFF1 Polyclonal Antibody



Catalog Number: 163552

---

Immunohistochemistry - TFF1 Polyclonal Antibody human stomach using TFF1 antibody at dilution of 1:200 (40x lens).  
Immunohistochemistry of paraffin-embedded