

MT-ATP8 Polyclonal Antibody

Catalog Number: 163554

Product name

MT-ATP8 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MT-ATP8

Preparation

Antigen: Recombinant protein of human MT-ATP8.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - MT-ATP8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MT-ATP8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 2min.



Immunofluorescence - MT-ATP8 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MT-ATP8 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

MT-ATP8 Polyclonal Antibody



Catalog Number: 163554

Immunofluorescence - MT-ATP8 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MT-ATP8 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.