

RNASE3 Polyclonal Antibody



Catalog Number: 163601

Product name

RNASE3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to RNASE3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-160 of human RNASE3 (NP_002926.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - RNASE3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using RNASE3 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - RNASE3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using RNASE3 antibody at dilution of 1:100 (40x lens).