

TIRAP Polyclonal Antibody



Catalog Number: 163630

Product name

TIRAP Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TIRAP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human TIRAP (NP_001034750.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - TIRAP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TIRAP antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - TIRAP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TIRAP antibody at dilution of 1:200 (40x lens).