

IPO5 Polyclonal Antibody



Catalog Number: 163646

Product name

IPO5 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to IPO5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human IPO5 (NP_002262.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - IPO5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IPO5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 3s.

Immunofluorescence - IPO5 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using IPO5 antibody . Blue: DAPI for nuclear staining.