

# FGR Polyclonal Antibody



Catalog Number: 163684

## Product name

FGR Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to FGR

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGR (NP\_005239.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



## Western blot - FGR Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FGR antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.