

LKB1 Polyclonal Antibody



Catalog Number: 163698

Product name

LKB1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LKB1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 134-433 of human LKB1 (NP_000446.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - LKB1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LKB1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry - LKB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using LKB1 antibody at dilution of 1:100 (40x lens).




LKB1 Polyclonal Antibody



Catalog Number: 163698


Immunohistochemistry - LKB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using LKB1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - LKB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human uterine cancer using LKB1 antibody at



dilution of 1:100 (40x lens).

Immunohistochemistry - LKB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophageal cancer using LKB1 antibody at dilution of 1:200 (40x lens).