

# THY1 Polyclonal Antibody



Catalog Number: 163700

## Product name

THY1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to THY1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human THY1 (NP\_006279.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:1000 - 1:2000


IF 1:50 - 1:200

## Validations




Western blot - THY1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using THY1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunofluorescence - THY1 Polyclonal Antibody

Immunofluorescence analysis of Mouse brain cells using THY1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



# THY1 Polyclonal Antibody



Catalog Number: 163700

---

Immunofluorescence - THY1 Polyclonal Antibody

using THY1 antibody at dilution of 1:100. Blue:  
DAPI for nuclear staining.

Immunofluorescence analysis of Rat brain cells