

REG3G Polyclonal Antibody



Catalog Number: 163709

Product name

REG3G Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to REG3G

Preparation

Antigen: Recombinant protein of mouse REG3G

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations



Western blot - REG3G Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using REG3G antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 10s.

Immunofluorescence - REG3G Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using REG3G antibody . Blue: DAPI for nuclear staining.