

# REST Polyclonal Antibody



Catalog Number: 163735

---

## Product name

REST Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to REST

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 760-1060 of human REST (NP\_001180437.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - REST Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using REST antibody .