

# IDH1 Polyclonal Antibody



Catalog Number: 163742

## Product name

IDH1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to IDH1

## Preparation

Antigen: A synthetic peptide of human IDH1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF, FC

## Dilutions


WB 1:500 - 1:1000

IHC 1:50 - 1:100

IF 1:20 - 1:50

FC 1:20 - 1:50

## Validations



Western blot - IDH1 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using IDH1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.