

# PSMA2 Polyclonal Antibody



Catalog Number: 163808

## Product name

PSMA2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PSMA2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human PSMA2 (NP\_002778.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - PSMA2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PSMA2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - PSMA2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PSMA2 antibody . Blue: DAPI for nuclear staining.