

AMPH Polyclonal Antibody



Catalog Number: 163848

Product name

AMPH Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to AMPH

Preparation

Antigen: Recombinant protein of human AMPH

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:25 - 1:100

Validations



Western blot - AMPH Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using AMPH antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.