

# CLCA4 Polyclonal Antibody



Catalog Number: 163861

## Product name

CLCA4 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CLCA4

## Preparation

Antigen: Recombinant protein of human CLCA4

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:100 - 1:500

IHC 1:10 - 1:50

## Validations



Western blot - CLCA4 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using CLCA4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.