

# DOK3 Polyclonal Antibody



Catalog Number: 163868

## Product name

DOK3 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to DOK3

## Preparation

Antigen: Recombinant protein of human DOK3

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:10 - 1:50

## Validations



Western blot - DOK3 Polyclonal Antibody

Western blot analysis of extracts of human notum skin, using DOK3 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.