

GNA11 Polyclonal Antibody



Catalog Number: 163877

Product name

GNA11 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GNA11

Preparation

Antigen: Recombinant protein of human GNA11

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - GNA11 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GNA11 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.