

HSF2 Polyclonal Antibody



Catalog Number: 163882

Product name

HSF2 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to HSF2

Preparation

Antigen: Recombinant protein of human HSF2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - HSF2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HSF2 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - HSF2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using HSF2 antibody at dilution of 1:200 (40x



lens).

Immunofluorescence - HSF2 Polyclonal Antibody

HSF2 Polyclonal Antibody



Catalog Number: 163882

Immunofluorescence analysis of A549 cells using HSF2 antibody . Blue: DAPI for nuclear staining.