

# OCT4 Polyclonal Antibody



Catalog Number: 163892

## Product name

OCT4 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to OCT4

## Preparation

Antigen: Recombinant protein of human OCT4

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



Western blot - OCT4 Polyclonal Antibody

Western blot analysis of extracts of K-562 cells, using OCT4 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.