

JUND Polyclonal Antibody



Catalog Number: 163915

Product name

Rabbit IgG

JUND Polyclonal Antibody

Applications

Specificity

WB, IHC

Human, Mouse, Rat

Dilutions

Antibody description

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Polyclonal antibody to JUND

Preparation

Validations

Antigen: A synthetic peptide of human JUND



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot - JUND Polyclonal Antibody

Clonality

Western blot analysis of extracts of Jurkat cells, using JUND antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.

Polyclonal

Ig Type