

PTGER2 Polyclonal Antibody



Catalog Number: 163936

Product name

PTGER2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PTGER2

Preparation

Antigen: A synthetic peptide of human PTGER2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:100 - 1:200

Validations




Western blot - PTGER2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PTGER2 antibody at 1:200 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - PTGER2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat lung using PTGER2 Antibody at dilution of 1:100



(40x lens).


PTGER2 Polyclonal Antibody



Catalog Number: 163936

Immunohistochemistry - PTGER2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using PTGER2 Antibody at dilution of 1:100 (40x lens).



Immunofluorescence - PTGER2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using PTGER2 antibody .