

# FOXN2 Polyclonal Antibody



Catalog Number: 163943

---

## Product name

FOXN2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to FOXN2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP\_002149.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - FOXN2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using FOXN2 antibody .